



Art.-No.: AO-A-12 (Single use)

Moorfield Forceps

Blunt tipped forceps for manipulating the conjunctiva.



Art.-No.: AO-C-22 (Single use) 105 mm

Colibri Style Toothed Forceps

1 into 2 teeth, 0,12mm with tying platforms



Art.-No.: AO-C-19 (Single use) 85 mm

Small Handle Straight Tying Forceps

Straight shafts with 4,5mm tying platforms



Art.-No.: AO-C-30 (Single use) 85 mm

Short Handle Capsularhexis Forceps

45° angled straight shafts. Tip to angled length 14,0mm
 Used to perform continuous tear circular or oval capsulatomies



Art.-No.: AO-C-20 (Single use) 85 mm

Small Handle Colibri Tying Forceps

1 into 2 teeth, 0,12mm with tying platforms



Art.-No.: AO-C-31 (Single use) 105 mm

Capsularhexis Forceps Straight

45° angled straight shafts. Tip to angled length 14,0mm
 Used to perform continuous tear circular or oval capsulatomies



Art.-No.: AO-C-21 (Single use) 100 mm

Straight Toothed Forceps

1 into 2 teeth, with 5,0mm tying platforms



Art.-No.: AO-C-32 (Single use) 105 mm

Capsularhexis Forceps Curved

45° angled gently curved shafts. Tip to angled length 14,0mm
 Used to perform continuous tear circular or oval capsulatomies



Art.-No.: AO-C-33 (Single use)

X Action Capsularhexis Forceps

45° angled gently curved shafts. Tip to angled length 12,0mm
Used to perform continuous tear circular or oval capsulatomies



Art.-No.: AO-C-51 (Single use)

Lens Folder

Folding 5-6mm sized optics of intra-ocular lenses (IOL).
Paddle style jaws with handle stop to prevent over-folding.



Art.-No.: AO-C-34 (Single use)

Round Bodied Capsularhexis Forceps

45° angled straight shafts. Tip to angled length 14,0mm
Used to perform continuous tear circular or oval capsulatomies



Art.-No.: AO-C-52 (Single use)

Lens Implantation Forceps

To implant the folded intra-ocular implant into the capsular bag.



Art.-No.: AO-C-49 (Single use) 115 mm

Straight Notched Forceps

0,25mm notched, with tying platforms



Art.-No.: AO-C-53 (Single use)

Locking Lens Introducer

Standard cross action with lock



Art.-No.: AO-C-50 (Single use) 115 mm

Curved Notched Forceps

0,25mm notched, with tying platforms



Art.-No.: AO-C-54 (Single use) 85 mm

St. Martins Forceps 1 in 2 teeth

1 into 2 teeth, with 5,0mm tying platforms



Art.-No.: AO-C-55 (Single use) 85 mm

Hoskins Forceps

0,25 mm notched, with tying platforms



Art.-No.: AO-G-02 (Single use) 12,5 cm

Mosquito Type Forceps

12,5 cm general holding forceps



Art.-No.: AO-C-71 (Single use) 105 mm

Straight Tying Forceps

Straight shafts with 4,5 mm tying platforms



Art.-No.: AO-A-11 (Single use)

Westcott Scissors

Curved scissors



Art.-No.: AO-C-72 (Single use) 100 mm

McPherson Tying Forceps

Angled 45° shafts with 4,5 mm tying platforms for tying knots and also manipulation of intra-ocular lens



Art.-No.: AO-C-80 (Single use)

Ong Type Scissors

Sharp pointed tips, angled blades



Art.-No.: AO-C-73 (Single use)

Kelmans Forceps

Angled shafts with tying platforms



Art.-No.: AO-C-81 (Single use)

Vannas Scissors

Straight scissors



Art.-No.: AO-G-01 (Single use) 95 mm

General Scissors Sharp Pointed

Sharp pointed 9,5 cm scissors



Art.-No.: AO-G-06 (Single use) 90 mm

General Scissors Blunt End


Blunt ended 9 cm scissors



Art.-No.: AO-G-14 (Single use) 130 mm

General Scissors Long Blunt End

Blunt ended 13 cm scissors



Art.-No.: AO-C-10 (Single use)

Kratz Type Paediatric closed wire speculum

10 mm blades



Art.-No.: AO-C-11 (Single use) 52 mm

Kratz Barraquer Wire Speculum

Speculum of steel wire with 13,5mm open blades

VE: 100 St. / pc.



Art.-No.: AO-C-12 (Single use)

Temporal Speculum

Speculum of steel wire, 15 mm blades



Art.-No.: AO-C-13 (Single use)

Closed Loop Kratz Speculum

Closed loop speculum of steel wire

VE: 100 St. / pc.



Art.-No.: AO-C-16 (Single use)

Pierce Type Speculum

Slide adjustable speculum



Art.-No.: AO-C-14 (Single use)

Sprung Speculum Paediatric

6 mm blades



Art.-No.: AO-C-15 (Single use)

Sprung Speculum Neonate

4 mm blade



Art.-No.: AO-R-10 (Single use)

Bladed Adjustable Speculum

Screw adjustable speculum with solid blades



Art.-No.: AO-R-11 (Single use)

Adjustable Speculum

Screw adjustable speculum, wire type



Art.-No.: AO-C-41 (Single use)

Y Nucleus Rotator

forked tip, blunt ends



Art.-No.: AO-C-42 (Single use)

Claw Phaco Chopper

curved chopper with two faceted chopping surfaces



Art.-No.: AO-C-43 (Single use)

Phaco Chopper

Ambidextrous nucleus chopper



Art.-No.: AO-C-44 (Single use)

Koch Type Chopper

Finger type nucleus chopper



Art.-No.: AO-C-45 (Single use)

Green Chopper
Angled nucleus chopper



Art.-No.: AO-G-15 (Single use)

Kuglen Hook
Iris hook and lens manipulator



Art.-No.: AO-C-46 (Single use)

Sinskey Lens Dialer / Sinskey Hook
For manipulating the intra-ocular lens.



Art.-No.: AO-R-21 (Single use)

Lasik Corneal Marker
Two blunt fins angled downwards and separated at 45° angle



Art.-No.: AO-C-47 (Single use)

Drysdale Nucleus Spatula
Paddle shaped tip



Art.-No.: AO-R-31 (Single use)

Flap Elevator
Blunt tipped 1mm spatula for manipulating the lasik corneal flap



Art.-No.: AO-C-48 (Single use)

Mushroom Nucleus Manipulator
for iris retraction



Art.-No.: AO-R-32 (Single use)

Universal Flap Lifter
Bevelled tip suitable for manipulating femto lasik flaps



Art.-No.: AO-R-34 (Single use)

Intra-Lasik Flap Separator

round, sharp tipped separator



Art.-No.: AO-R-40 (Single use)

CK Lasik Marker (target center)

Traget center to facelitate accurate marking of the cornea



Art.-No.: AO-R-41 (Single use)

Lasek Alcohol Well 9 mm diameter

Liquid holding well for alcohol application during surgery



Art.-No.: AO-R-42 (Single use)

Lasek Alcohol Well 8 mm diameter

Liquid holding well for alcohol application during surgery



Art.-No.: AO-R-43 (Single use)

Lasek Alcohol Well 7 mm diameter

Liquid holding well for alcohol application during surgery



Art.-No.: AO-R-51 (Single use)

Epi-peeler

For manipulating the epithelium in lasek procedures



Art.-No.: AO-R-60 (Single use)

Hoe Type Epi-peeler (D Shaped)

For manipulating the epithelium in lasek procedures



Art.-No.: AO-R-61 (Single use)

Hoe Type Epi-peeler (Square Shaped)

For manipulating the epithelium in lasek procedures



Art.-No.: AO-R-71 (Single use)

Lasik Flap Forceps

For manipulating the flap in lasik procedures



Art.-No.: AO-CR-01 (Single use)

Curette Small

Chalazion curette - 1 mm diameter



Art.-No.: AO-R-90 (Single use)

Lasik Cannula (single hole at tip)



Art.-No.: AO-CR-02 (Single use)

Curette Medium

Chalazion curette - 1,5 mm diameter



Art.-No.: AO-CS-01 (Single use)

Chalazion Clamp small

11mm diameter, round plate



Art.-No.: AO-CR-03 (Single use)

Curette Large

Chalazion curette - 2 mm diameter



Art.-No.: AO-CS-02 (Single use)

Chalazion Clamp medium

Oval plate, 23 x 16mm



Art.-No.: AO-PD-01 (Single use)

Punctal Dilator Sharp



Art.-No.: AO-PD-02 (Single use)

Punctal Dilator Blunt



Art.-No.: AO-A-12 (Single use)

Moorfields Forceps

Blunt tipped forceps for manipulating the conjunctiva



Art.-No.: AO-G-04 (Single use)

Braunstein Fixed Caliper

3,5 / 4 mm



Art.-No.: AO-G-05 (Single use)

Castroviejo Adjustable Caliper

Caliper measures from 0mm to 20mm in 0,5mm increments with scale reading on one side



Art.-No.: AO-G-02 (Single use)

Mosquito Type Forceps 12,5 cm

12,5 cm general holding forceps



Art.-No.: AO-G-08 (Single use)

Jaffe Type Lid Retractor Top LID

14 mm blades (x2)



Art.-No.: AO-G-03 (Single use)

Trabeculectomy Punch

Bullet shaped tip for easy insertion into sclero-corneal tunnel



Art.-No.: AO-G-09 (Single use)

Jaffe Type Lid Retractor Bottom LID

9 mm blades (x2)



Art.-No.: AO-G-10 (Single use)

Jameson Strabismus Hook

5 mm hook



Art.-No.: AO-G-15 (Single use)

Kuglen Hook

Iris hook and lens manipulator



Art.-No.: AO-G-11 (Single use)

Jameson Hook Extended

8 mm hook



Art.-No.: AO-C-23 (Single use)

Neonate Scleral Depressor Bulbous end



Art.-No.: AO-G-12 (Single use)

Fine Thornton Fixation Ring small

Swivel head 13mm/11mm cut out



Art.-No.: AO-C-24 (Single use)

Neonate Scleral Depressor Hammer shape



Art.-No.: AO-G-13 (Single use)

Fine Thornton Fixation Ring standard

Swivel head 16mm/14mm cut out



Art.-No.: AO-C-61 (Single use)

Needle Holder

Curved, 6 mm delicate jaw

Argonoptics GmbH & Co KG
Zum Stadtgraben 7
45721 Haltern am See

www.argonoptics.de

Tel. 02364 50 88 70